

Single Embryonic Mouse Primary Cortical Neuron Culture--Dissociated via Trypsin

Day 1

Poly-D-Lysine coating:

- 1.) Add 333.35mL of sterile ddH₂O to 50mg of PDL
- 2.) Mix and filter thru 0.22 µm filter.
- 3.) Aliquot and store @ -20°C
 - a.) Use within 1-2 weeks
 - b.) Only freeze / thaw once...

To coat plates:

- 1). Thaw PDL at 37°C
- 2). Filter thru 0.22 µm filter
- 3). Plate and incubate at 37°C overnight

Plating volumes:

96 well plates	50-60 µL/well.
4 well slides	300 µL/well
24 well plates	500 µL/well
6 well plates	2 mL/well
T25 flasks	4 mL/flask

Day 2

Wash / Dry Poly-D-Lysine Coated Plates

- 1.) Vacuum off PDL and wash twice with sterile ddH₂O
- 2.) After the second wash is removed, allow the plates to dry completely (uncovered in the hood)

Prepare Mouse Laminin Substrate:

- 1.) Add 40 mL of base media to 1.1 mg of mouse laminin
* Do not filter *
- 2.) Coat PDL plates with laminin (same volumes as for PDL) and place in incubator overnight

Day 3

Wash Laminated plates

- 1.) Vacuum off laminin and wash twice with an equal volume of base media.
- 2.) After the second wash is removed, replace a third time with base media and place in incubator until time of plating.

Dissection / Preparation:

- 1.) Warm dissection media (HBSS w/o calcium or magnesium) and place in appropriate dishes.
- 2.) Warm plating media in 37°C H₂O bath.
- 3.) Sacrifice the mother by cervical dislocation.
- 4.) Spray lower abdomen with 70% EtOH and cut through skin and muscle with a pair of scissors
exposing the uterus, intestines and embryos.
- 5.) Cut embryos from amniotic sac.
- 6.) Remove brain from skull.
- 7.) Dissect the cortex from each lobe.
- 8.) Using a spring scissors, cut each cortex into about 3 pieces
- 9.) Pick up the cortical pieces along with 1 mL of the HBSS using a blue tip and place into a 15 mL conical tube.

- 10.) Add 100 μ L of 0.25% Trypsin to the 1 mL of dissection medium
Final trypsin concentration = 0.025%
- 11.) Place tissue / enzyme in 37°C incubator for 30 minutes.
- 12.) After the 30 minute digestion is complete, spin the tissue at 1500 rpm for 3 minutes.
- 13.) Remove the supernatant and re-suspend in 1 mL of plating media.
- 14.) Count cells on hemacytometer and dilute to 0.5×10^6 cells/mL

Feeding / Inhibition:

- On DIV5, remove 50% of the media and replace with fresh plating media along with Ara-C at a final concentration of 1 μ M.
 - Aliquots of 10 mM AraC in ddH₂O are stored @ -20°C
 - 20 μ l of 10 mM AraC is added to 100mL of feeding media = 2 μ M
 - 1/2 change of media = final of **1 μ M**.
 - Note: Never do a 100% media change.
- After the DIV5 inhibitory feeding, feed the cells every 3 days (without Ara-C) using the 50% exchange method and replacing with feeding media (plating media but without serum, 2% B27 is supplemented instead).

Solutions:

Neurobasal Plating Media

- 1.) Defrost cFBS and HS in 37°C H₂O bath.
- 2.) From 500 mL Neurobasal bottle, remove 50 mL.
- 3.) Add 25 mL of cFBS and 25 mL of HS
- 4.) Filter through 0.22 μ m filter system, label and store @ 4°C
 - * When ready to plate / feed, add FRESH L-glutamine (200 mM stock, 2mM final) and Pen-Strep (100% stock, 0.5 to 1% final).

NB/B27 Feeding Media

- 1.) Defrost B27 in 37°C H₂O bath.
 - 2.) From 500 mL Neurobasal bottle, remove 10 mL.
 - 3.) Add 10 mL of B27
 - 4.) Filter through 0.22 μ m filter system, label and store @ 4°C.
- When ready to feed, add FRESH L-glutamine (200 mM stock, 2mM final) and Pen-Strep (100% stock, 0.5 to 1% final).

Product Ordering Information:

0.25% Trypsin, 15050-065, Invitrogen
 Ara-C, C6645, Sigma
 B27 Supplement, 17504044, Invitrogen
 FBS (Dialyzed), SH30071.03HI, HyClone
 L-glutamine, 25030-081, Invitrogen
 Neurobasal Media, 21103-049, Invitrogen
 NUNC 6 well plates, 140675, Fisher
 Costar 24 well plates, 3526, Fisher
 Poly-D-Lysine, P0899, Sigma