

Hackam Lab

Whole Mount Protocol for Intestine, Lung, Liver and Brain (Embryo and postnatal tissues)

*****Keep plate wrapped with foil through all steps if tissues are fluorescent labeled*****

1. Harvest tissues, mark outer surface with INIDA INK before fixing in 4% Paraformaldehyde (PFA)
(For intestine—keep the intestine on paper towel, add a drop of INDIA INK, wait for couple of seconds, remove/dry excess INDIA INK, cut intestine longitudinally to open the lumen).
2. Place tissues in wells (24-well plate) containing ~1mL 4%PFA, and fix at 4°C (on rotator) for 2 hours.
3. Remove PFA, refill each well with 1 ml 1x PBS, allow washing on rotator for 10 minutes at room temperature.
4. Remove PBS and refill each well with 1 ml 1x PBS, allow washing on rotator for 10 minutes at room temperature.
5. Remove PBS and refill each well with 1ml 1x PBS (**exact volume**), allow washing on rotator for 10 minutes at room temperature.
6. Remove 250 μ L of PBS and add 250 μ L methanol back to each well, allow washing on rotator for 15 minutes at room temperature *(wells now contains 750 μ L PBS: 250 μ L Methanol).*
7. Remove 500 μ L PBS: methanol and add 500 μ L methanol to each well, allow washing on rotator for 15 minutes at room temperature *(wells now contains 500 μ L PBS: 500 μ L methanol).*
8. Remove 750 μ L of PBS: methanol solution and add 750 μ L methanol to each well, allow washing on rotator for 15 minutes at room temperature *(wells now contains 250 μ L PBS: 750 μ L methanol).*
9. Remove whole 1 ml of PBS: methanol solution and add 1 ml 100% methanol to each well, allow washing on rotator for 15 minutes at room temperature *(wells now contains 100% methanol).*
(You may STOP and store for later; just seal with plastic wrap to avoid methanol evaporation, cover with aluminum foil and store in -20°C freezer).

10. Remove methanol and refill each well with 1 ml 1x PBS, allow washing on rotator for 15 minutes at room temperature.
11. Remove PBS and refill each well with 1 ml 1x PBS, allow washing on rotator for 15 minutes at room temperature.
12. Remove PBS and refill each well with 1 ml 1x PBS, allow washing on rotator for 15 minutes at room temperature.
13. **BLOCKING:** Block with 10% donkey serum in PBST for 1 hour at room temperature.
14. **PRIMARY ANTIBODIES:** Add primary antibodies (you could use multiple antibodies) (1:250 dilution in 1% donkey serum/PBST) incubate in 4°C refrigerator on rotator 1–2 nights (*2 nights sometime better*).
15. Wash with PBST on rotator in 4°C for 4–5 hours changing PBST **at least 5 times**.
16. **SECODARY ANTIBODIES:** Add secondary antibodies (1:1000 dilution in 1% donkey serum/PBST, incubate in 4°C refrigerator on rotator overnight (*24 hours or more*)).
17. Wash with PBST 4°C refrigerator on rotator overnight for 4–5 hours changing PBST **at least 5 times**.
18. Add nuclear stain (DAPI for 10–15 minutes) directly to secondary if needed solution; incubate for 30 minutes at room temperature.
19. Wash 6 times with PBST at 4°C each time incubate for 10 minutes on rotator.
20. Mount with Gelvatol mount media (*use dissection microscope to make sure lumen of intestine is facing up toward coverslip; tighten the coverslip and slide using small paper clips*).

Reagents:

- 4% PFA—Dilute one capsule of 16% PFA (Electron microscopy Sciences Cat # 15710) in 1x PBS.
- PBST—0.1% Triton-X 100 in 1xPBS.
- 10% donkey serum in PBST—add 100µL donkey serum in 900µL PBST.
- DAPI solution: DAPI (4',6-Diamidino-2-Phenylindole, Dilactate) (BioLegend cat # 422801)—To make 5 mg/ml (10.9 mM) solution, dissolve entire contents of vial in 2 ml of deionized water. Store between 2°C and 8°C, protected from light. For long-term storage, aliquot and store at -20°C.