



## HOME ECLIPSE DEVICE

Some medications are stored at room temperature while others are refrigerated. If your medication is refrigerated, be sure to remove your medication from the refrigerator, as instructed. If you forget to take your scheduled dose out of the refrigerator in advance, it will not hurt you to administer it cold. However, it may feel cold as it is infusing and it may slow down the rate of infusion.

1. Wash Hands or Hand Sanitize
2. Prepare a clean work surface area.
3. Gather your supplies: alcohol wipes, 2 saline syringes (white top), 1 Heparin syringe\_\_\_\_\_ (if applicable) and IV Eclipse medication. Compare the medication label to the prescription in your home folder. Make sure it has your name, the correct medication, correct dose and the “do not use after” date is not expired.
4. Re-clean your hands.
5. Unwrap the flushes and place on clean surface area.
6. Use SASH method: **S** – Saline (white cap)  
**A** – Administer  
**S** – Saline (white cap)  
**H** – Heparin (\_\_\_\_\_ cap)
7. Clean the end cap or hub with an alcohol wipe for 15 seconds (“scrub the hub”).
8. Attach the saline syringe; open the catheter clamp and flush.
9. Discard the empty syringe.
10. Clean the end cap or hub again with a new alcohol wipe for 15 seconds (“scrub the hub”).
11. Remove the protective cover from the end of the tubing on your Eclipse ball and attach to the end of your catheter. Open the clamp on the Eclipse ball.
12. Allow the medication to infuse for the prescribed time. When it is finished, it will have no bulges along the inner column.
13. When the infusion is complete, close the clamp on the Eclipse ball and disconnect the tubing from your catheter, leaving the end cap in place.
- 14.

15. Clean the end cap or hub again with a new alcohol wipe for 15 seconds (“scrub the hub”).
16. Flush your catheter with saline, as you did above.
17. Clean the end cap or hub again with a new alcohol wipe for 15 seconds (“scrub the hub”).
18. Flush your catheter with Heparin (if applicable) and then clamp your catheter.
19. Discard used syringes and the Eclipse ball into the garbage, or as instructed.
20. Re-clean your hands.