

## **IVIG TAKE DOWN**

When your infusion is complete the pump will either alarm END PROGRAM (the pump has administered the total volume amount) or AIR/UPSTREAM OCCLUSION (all the medication in the bag has infused but the pump thinks there is still more volume left). For either alarm, open the pouch and check to see if the bag is empty. If there is any remaining fluid in the bag please call us for further instruction (day 410-288-8105 #1 and weekends/afterhours 410-288-8000 #4).

If the bag is empty, your infusion is complete.

1. Press the RED STOP BUTTON one time. This will stop your pump.
2. Press and hold down the GREY ON/OFF BUTTON until the back bar scrolls across the screen. This will turn your pump off.
3. Wash hands/Hand sanitize.
4. Disconnect the tubing from your IV catheter, leaving needless end cap/hub on the IV line.
5. If you are scheduled for another treatment tomorrow you will need to flush your peripheral IV with Saline and Heparin. If you are not scheduled for another treatment tomorrow, it is OK to remove the peripheral IV at this time.
6. Clean the end cap or hub with an alcohol wipe for 15 seconds (“scrub the hub”)
7. Flush your catheter with 5 ml of Normal Saline.
8. Clean the end cap or hub with an alcohol wipe for 15 seconds (“scrub the hub”)
9. Flush your catheter with 5 ml of heparin and then clamp your catheter.
10. Discard used syringes, medication, and tubing in the garbage.
11. Your end cap/hub does not need a cover.